




Lab Facilities

- AwiLAB Gasanalysis
- AwiLAB Digester






AWILAB GASANALYSIS



AWILAB GAS ANALYSIS

Our product family AwilAB was especially developed for biogas application but can also be used for various researches in the gas sector.

The AwilAB Gasanalysis guarantees smooth operating processes and with that your research. In combination with our AwifLEX system, the gas production and composition can be controlled and evaluated during the substrate fermentation process.



SCOPE OF APPLICATION


Batch/ discontinuous test

- gas production less than 4 l/ hr
- determination of
 - daily gas production
 - start/ end of fermentation process
 - time of maximum gas production
 - ...
- comparison of different substrates with different TS contents or mixing ratios
- comparison from process changes during test (temperature, pressure, etc.)
- effect of the addition of additives (trace elements, etc.)
- effect of gas volume on organic load
- detection of
 - process disturbances in long-term operation
 - changes in gas composition

➔ To answer scientific questions, a reference is required in addition to the experimental digester.


SET-UP

AwiLAB 1




SET-UP


Gasbag 1 Gasbag 2 Gasbag 3 Gasbag 4




SET-UP



SET-UP




SET-UP



CABLE WITH 3.5 STEREO JACK ENDINGS



FEATURES

- stainless steel housing (WxHxD = 500x500x210 mm) for wall mounting completely assembled.
- connections for max. 21 Milligascounter (MGC) for triple determination via 3.5 mm jack plug (3-pin) or
- connection of max. 7 drum-type gas meters via 3.5 mm jack plug (3-pin) or
- combination of MGC and drum-type gas meters
- max. 14 process connection (threaded hose connections made of stainless steel 4/6 mm)
 - 7 connections for measuring point
 - 7 connection for gas bags
- control of encapsulated analysis valves (changeover to exhaust air with full bag)
- max. 7 LED indicator lights for status display
- 1 measuring point extension to gas analysis system for pressure compensation (N2)
- data connection to Awite process analysis system via RS485
- AwiCore electronic module consisting of:
 - controller module
 - AWICORE CBase fitted with:
 - 2x 8x digital output plug-in card
 - 3x 8x digital input plug-in card
- maximum combination with
 - max 63 pcs. Milligascounters for triple determination
 - max 21 pcs. Drum-type gas meters or
 - combination of MGC and drum-type gas meters



FEATURES

- recording the gas volumes produced per fermenter
- recording the current filling levels of the gas bags per fermenter
- automatic release for measurement at a certain adjustable volume
- automatic emptying of the bag after measurement via pressure sensor
- reset function for individual counter and filling level for bag
- input of counter values (milligascounter/ drum-type gas meter)
- input of calibration factors
- visualization on 7" touch panel of the gas analysis system Awiflex
- with milligascounter : Fault detection when the MGC is stuck



VISUALISATION & MENU NAVIGATION

Screenshot of a user interface showing a navigation menu and various settings.

The top navigation bar includes:

- State
- Recent Values
- History
- Admin
- Help
- Contact
- Settings
- Measurement Plan
- Cooler
- Cooler 2
- Service

On the left side, there is a vertical sidebar with the following information:

- 2372
- 16:00

The main content area contains several buttons and a slider:

- Start Measuring
- Controller Settings
- Calibration
- Output Report
- Change Userlevel (with a value of 100)
- Userlevel
- Expert Settings
- AwiLAB (highlighted with a yellow circle)
- Power Off
- USB-Stick
- Language

A yellow arrow points from a callout box at the bottom to the "AwiLAB" button.

Callout box text: special settings for AwiLAB


VISUALISATION & MENU NAVIGATION




VISUALISATION & MENU NAVIGATION

Batch

- triple- determination
- 3 reactors (A/B/C) fill one gasbag




VISUALISATION & MENU NAVIGATION




VISUALISATION & MENU NAVIGATION

Continuous stirred digesters (one gasmeter)




EVALUATION

Application: Batch




EVALUATION

Application: Batch



EVALUATION

CH₄, CO₂ and O₂ values



AWILAB DIGESTER



AWILAB DIGESTER

Our Awilab Digester is an all-in lab facility for your research.

„Stand Alone“



„Modular“



DIGESTER
MODULE

BASIC
MODULE

EQUIPMENT

The AwiLAB Digester ...

- stand-alone or modular system - up to 7 digester modules
- compact & space-saving
- variable digester sizes/ stainless steel digester especially for biogas application/ completely insulated
- sight glasses for visual process monitoring
- process-optimised stirrer technology/ high-performance stirrer/ simmering seal
- customized stirrer configuration/ specific interval stirring
- removable stirrer for safe and easy removal
- individual heating program
- continuous monitoring and adjustment of temperature
- sampling at the most homogeneous zone
- special seals on the entire fermenter
- integrated self-pressure test for leak-test of the digester
- integrated overpressure protection
- operation und visualisation with 7-TFT-Touch-Display


EQUIPMENT

The AwiLAB Digester ...

- fully automatic phase monitoring
- automatic measuring point switchover
- fully automatic measurement of gas quality
- automatic normalisation of the gas volume
- gas storage bag (5l); pressure equalisation fermenter (2 l) and after measurement (15 l)
- time-controlled emptying of the gas bags
- continuous monitoring of trials - also via remote access! (AwiMobileControl)
- online process measuring and control technology/ various ports for measurement technology
- central data collection and storage
- graphical and tabular analysis (AwiCharts/ AwiControl/ AwiView)
- extensive alarm and error messages
- no additional software licenses

SET-UP


PROCESS
OPTIMISATION



USER-FRIENDLY

INDIVIDUALLY
CONTROLLABLE


COMPACT



HIGH
QUALITY
STANDARD

1. Digester
2. Control cabinet incl. display
3. Electrical cabinet gas collector
4. Drum-type gas meter
5. Overpressure protection
6. Gas analysis system
7. Gas compensation bag
8. Inert gas bag
9. Gas storage bag analysis gas


SET-UP




- | | |
|----------------------------|-------------------------|
| 1. Digester | 9. Optional sensor |
| 2. Filling mechanism | 10. Optional sensor |
| 3. Overpressure protection | 11. Optional connection |
| 4. Condensate separator | 12. Optional connection |
| 5. Sampling port | |
| 6. Emptying port | |
| 7. Engine with stirrer | |
| 8. Temperature sensor | |

optional further installations possible

SET-UP



SET-UP



VISUALISATION & MENU NAVIGATION

The screenshot illustrates the Awite software interface for monitoring a gas analysis digester. The interface is divided into several sections:

- Header:** State, Settings, Gasanalyzer, History, Contact, Admin, Tab, Tab.
- Left Sidebar:** Shows a device icon, the ID 1880, and a U-Level button.
- Time and Date:** 08:58
- Central Area:**
 - Ambient Sensors:** Temperature / °C: 21.3, Abs. pressure / mbar: 949.
 - Digester Sensors:** Temperature / °C: 22.7, Pressure / mbar: -0.3.
 - Phase:** Startup, Outgassing, Operation.
 - Show more information:** A blue button at the bottom of the sidebar.
- 3D Model:** A detailed 3D rendering of the gas analysis digester unit.
- Gas Composition Data:**

Gasanalysis digester	
Methane (CH ₄) / %	0
Oxygen (O ₂) / %	21.65
Carbon dioxide (CO ₂) / %	0.1
Hydrogen sulfide (H ₂ S) / ppm	3
Hydrogen (H ₂) / ppm	0
- Gas Counter Data:**


Gascounter normalized according to ambient sensors 0°C/1013mbar	
Flow / l/h	0.000
Volume / l	110.0
F per m ³ substrate / m ³ /h	0.000
Water content / %	2.7

Annotations with orange circles and arrows highlight specific features:

- AMBIENT CONDITIONS:** Points to the Ambient Sensors section in the sidebar.
- INDIVIDUAL CONFIGURED:** Points to the Gas Composition data table.
- GAS COMPOSITION:** Points to the Gas Composition data table.
- STANDARDISATION:** Points to the Gas Counter data table.

VISUALISATION & MENU NAVIGATION

MORE DETAILS



Other		Gascounter		Gasanalysis	
Stirrer		Normalized according to ambient sensors 0°C/1013mbar		Digester	
Setpoint / %	56	Flow / l/h	0.000	Methane (CH4) / %	0
Current / A	0.54	Volume / l	14.405	Oxygen (O2) / %	21.6
Ambient sensors		F per m³ substrate / m³/h	0.000	Carbon dioxide (CO2) / %	0.6
Temperature / °C	26.0	Water content / %	3.4	Hydrogen sulfide (H2S) / ppm	1
Abs. pressure / mbar	988	Not normalized Raw data		Hydrogen (H2) / ppm	0
Digester		Flow / l/h	0.000	Ambient air	
Temperature / °C	40.9	Volume / l	16.740	Methane (CH4) / %	0
Pressure / mbar	-0.1	F per m³ substrate m³/h	0.000	Oxygen (O2) / %	21.63
Collecting/Measuring gas		Gasbag (max. 4.2 l)		Carbon dioxide (CO2) / %	0
Inactive		Volume / l	0.0	Hydrogen sulfide (H2S) / ppm	3
Phase				Hydrogen (H2) / ppm	0
Startup				Actuators	
Outgassing				Stirrer	
Operation				Heater	
Show less information					

VISUALISATION & MENU NAVIGATION

INDIVIDUALLY
CONTROLLABLE

	State	Settings	Gasanalyser	History	Contact	Admin	Tab	Tab
	Experiment	Stirrer	Gascounter	Other				
1880								
	Experiment				Startup		Outgassing	Operation
	Start	Stop			Use	Skip	Use	Skip
	Stirring	Interval	On	Off				
	Interval runtime / min				1		0.2	1
	Interval pausetime / min				0		0	0
	Setpoint (0..100) / %				56		56	56
	Inactivity alarm after / min				5		1	1
	Heating	Auto	On	Off				
	Set temperature / °C				22.6		41.0	50.0
	Hysteresis / °C				0.2		0.5	0.2
	Increase by / °C				0.0			
	Increase within / h				0.0			
	Max. deviation alarm +- / °C				0.0		0.0	0.5
	Activity alarm after / min				1440		1440	1440
09:20	Other							
	Duration / h m s				0 00 00		0.0	3 21 43
	Countdown / h m s				-49 -55 -25		0 00 00	
	Analyze Gas?				Yes	No	Yes	No

VISUALISATION & MENU NAVIGATION

**INTEGRATED
PRESSURE TEST**

State **Settings** **Gasanalyser** **History** **Contact** **Admin**

Experiment **Stirrer** **Gascounter** **Other**

Pressure test

Start test	Manual	Off
Actual pressure / mbar		0.1
Set pressure / mbar		18.0
Allowed loss / mbar		2
Equalisation time / s		35
Maximum attempts		3
Holding time / s		130
Max. runtime pump / s		700
Last test: SUCCESS		
Start pressure / mbar		0
Equalisation countdown / s		0
Test countdown / s		0

Other


Digester Volume / l	50.0
---------------------	------

1880
U-Level

14:58


EVALUATION EXAMPLE AWICHARTS

CH4 / CO2



EVALUATION EXAMPLE AWIVIEW

CH4 / CO2



SUMMARY



compact, modular

reliable in process
all-in system



extendable
individual

multi-functional
extendable
smart/ clever
extensive/ wide



custom-made
easy to use
minimum operating effort
maximum operating safety



remote monitored
completely controllable
fully monitored
no additional licenses necessary



individual
high quality equipment
visualised



REFERENCES



BILDUNGSZENTRUM TRIESDORF



AWITE BIOENERGIE GMBH
INFO@AWITE.DE | WWW.AWITE.COM
TEL. +49 (0)8761 / 72162-0

